

INDUSTRIAL / NON-SURGICAL

“FOLDABLE” N95 RESPIRATOR

Manufactured in USA with globally sourced material.

Part Number: PT-F01

Description: Foldable N95 rated respirator with moldable nose clip

Location: Manufactured in Ohio, USA

Brand: Patriot Mask



Made in the USA with globally sourced materials!

KEY FEATURES:

- ≥95% BFE (Bacterial Filtration Efficiency) according to ASTM F2101.
- ≥95% Filter Efficiency according to ICS Laboratories Inc.
- Lightweight construction promotes longer wear time.
- Moldable nose clip will help provide custom fit and secure seal.
- One size fits most.
- Latex free: Yes
- Model: Foldable.

DO NOT USE:

- For gases and Vapors, including oil aerosols, etc.
- If modified or used for its unintended purpose
 - > Expiration date of 3 years after manufacture date

MATERIAL COMPOSITION:

- **Outer Layer:** Spun bound cloth/non-woven fabric.
- **Filter:** Electrostatic hot air cotton.
- **Lining:** Melt blown non-woven fabric.
- **Inner Layer:** Spun bound cloth/non-woven fabric.
- **Nose Clip:** Aluminum moldable.
- **Nose Foam:** Spandex filament sponge bar.
- **Headband:** Nylon string.
- **Breathing valve:** N/A

DO USE FOR:

- Medical clinics and operating rooms
- Laboratory procedures
- Infection control requirements
- Emergency planning and preparedness

INDUSTRIAL / NON-SURGICAL

“FOLDABLE” N95 RESPIRATOR


Manufactured in USA with globally sourced material.

Part Number: Patriot F-01
Description: Foldable N95 rated respirator with moldable nose clip
Location: Manufactured in Ohio, USA
Brand: Patriot Mask

PACKAGE INFORMATION

Description	PCS/Box	Box/Carton	Pcs/Pkg	Carton/Pallet	Total Quan
Patriot F-01	40	20	800	24	19,200

DIMENSIONS



FITTING INSTRUCTIONS:

1. Hold the respirator in hand and with the nosepiece at your fingertips. Allow headbands to hang freely below hand.
2. Wear the respirator firmly against your face with the nosepiece on the bridge of your nose.
3. Stretch and pull the lower headband over the head and position below your ears. Stretch and pull the top headband on the back of your head above your ears.
4. Press soft metal to conform snugly around your nose.
5. Seal check.
 - To test the fit of a respirator without an exhalation valve, cup both hands over the respirator and exhale sharply.
 - To test the fit of a respirator with an exhalation valve, cup both hands over the mask and inhale sharply. A negative pressure should be detected inside the respirator.
 - If air flow is felt in the nose area, re-adjust/tighten the nose clip.
 - If air flow felt around the edges of the respirator, re-position the respirator / head strap to achieve a better fit.
6. Change the mask immediately if breathing becomes difficult or mask becomes damaged or distorted.
7. Change the respirator if a proper face seal cannot be achieved.
8. Careful observance of these instructions is an important step in safe respirator use.

WARNING:

95% efficient against particulate aerosols free of oil. This respirator helps protect against certain particulate contaminants but does not eliminate exposure to or the risk contracting disease or infection. Misuse may result in sickness or death.