

SPECIFICATION SHEET

Part Number:	PT-N95F-07
Description:	Foldable N95 rated respirator with moldable nose clip
Location:	Manufactured in Ohio, USA
Brand:	Patriot Mask



NIOSH
TC-84A-9282



Key Features:

- Meets NIOSH inhalation procedure testing TEB-APR-STP-0007.
- 99% BFE (Bacterial Filtration Efficiency) according to ASTM F2101.
- 99.7% Filter Efficiency according to ICS Laboratories Inc.
- Berry Amendment (10 U.S.C. 2533a) compliant.
- Lightweight construction promotes longer wear time.
- Moldable nose clip will help provide custom fit and secure seal.
- One size fits most.
- Latex free: Yes.
- Model: Foldable.

Material Composition:

- **Outer Layer:** Spun bound cloth/non-woven fabric
- **Filter:** Electrostatic hot air cotton
- **Lining:** Melt blown non-woven fabric
- **Inner Layer:** Spun bound cloth/non-woven fabric
- **Nose Clip:** Aluminum moldable
- **Nose Foam:** Spandex filament sponge bar
- **Headband:** Nylon string
- **Breathing valve:** N/A



Approvals and Standards:



TC-84A-9282



Manufactured in the USA

Do NOT Use:

- For gases and Vapors, including oil aerosols, etc.
- For repeated use following NIOSH guidelines
- If modified or used for its unintended purpose
- > Expiration date of 3 years after manufacture date

DO Use For:


- Medical clinics and operating rooms
- Laboratory procedures
- Infection control requirements
- Emergency planning and preparedness

Package Information:


Description	PCS/Box	Box/Carton	Pcs/Pkg	Carton/Pallet	Quantity
PT-N95F-07	40	20	800	32	25,600

Unit: mm

Front view



Side view



Back View

